

FOM-02M FiberOptoMeter II

The "optical multimeter" for probing brain circuits

NEW

FOM-02M/2 Dual Excitation/Dual Fluorescence Detection

FOM-02M Dual Excitation/Single Fluorescence Detection

Optogenetic Stimulation &
Fluorescence Measurement
via the Same Fiber



Features

- ⇒ **Dual-band excitation:**
450 – 490 nm (OGB1, GCaMP, ChR2)
550 – 590 nm (ArchT)
Fluorescence recording in 500 – 540 nm band.
- ⇒ **short pulses at high power** (optogenetics)
on top of
- ⇒ **constant low power illumination** (fluorescence)
- ⇒ **Ultrafast photodetector:**
high robustness against high-dose irradiation
detecting fluorescence signal immediately after
optogenetic stimulation
- ⇒ can also be used as **supplemental
photodetector** channel for **multiphoton
microscopy**.
- ⇒ compatible with multi-mode SMA fibers
incl. **non-magnetic fibers** for **fMRI**.

Applications


- ⇒ Localized measurement of physiological
fluorescence signals (e.g. Ca^{2+} signals) from
multiple sites in **freely behaving animals
including primates**
- ⇒ **Multi-color optogenetic** stimulation
- ⇒ **Near simultaneous optical measurement of
cellular activity and optogenetic stimulation**
- ⇒ **Chronic monitoring** of physiological signals in
disease model animals and drug testing
- ⇒ **Enhanced localized photodetection for
deep-tissue multiphoton imaging**
- ⇒ **Combined fluorescence recording with fMRI
imaging**



Examples of Fluorescence Measurement *via* a 200 μm Fiber

Upper trace:

Slow calcium waves (isoflurane 1.5%)
spontaneous activity (200 μm fiber)



Lower trace:

Same measurement as above,
visually evoked (*) and
spontaneous slow calcium waves



Ca²⁺-Traces

Ca²⁺ fluorescence indicator OGB-1 was injected into the visual cortex of a mouse.
Data kindly provided by Dr. A. Stroh and M. Schwalm.

REF: **Justus et al.** (2016), Glutamatergic synaptic integration of locomotion speed via septoentorhinal projections. *Nat. Neurosci.* <http://dx.doi.org/10.1038/nn.4447>


Fuhrmann et al. (2015), Locomotion, theta oscillations, and the speed-correlated firing of hippocampal neurons are controlled by a medial septal glutamatergic circuit. *Neuron*, <http://dx.doi.org/10.1016/j.neuron.2015.05.001>

Adelsberger et al. (2014), In vivo calcium recordings and channelrhodopsin-2 activation through an optical fiber. *Cold Spring Harb Protoc.* <http://dx.doi.org/10.1101/pdb.prot084145>

Stroh et al. (2013), Making Waves: Initiation and Propagation of Corticothalamic Ca²⁺ Waves In Vivo. *Neuron*, <http://dx.doi.org/10.1016/j.neuron.2013.01.031>

Technical Specifications:

Available central wavelengths (LED):	450 nm, 470 nm, 570 nm, 617 nm, 625 nm
Excitation band:	450 – 490 nm (blue) 540 – 590 nm (orange)
Max. light power at 200 μm fiber tip:	110 mW/mm ² (blue) 65 mW/mm ² (orange)
Fiber interface:	SMA or FC/PC
Fluorescence detection band:	500 – 540 nm (green)
Fluorescence signal output filter:	200 Hz low-pass (add-on block)
Fluorescence signal output range:	0 ~ -5 V (negative polarity)
Dimensions:	Requires 2 module slots (FOM-02M) or 5 module slots (FOM-02M/2) of EPMS chassis (npi electronic)
Power consumption:	20 W



L1 - L4: collimation lenses
F1 - F4: bandpass filter plates
DM1 - DM2 : dichroic mirrors
PMT: photomultiplier module

The optogenetic stimulation and fluorescence measurement via the same fiber was developed in the Konnerth Lab, TU München, Germany (Adelsberger et al *Nat Neurosci.* 2005, (8):988-90). The FOM-02M was designed in collaboration with Dr. Hongbo Jia, Suzhou Institute of Biomedical Engineering and Technology.

General:

npi electronic GmbH
Phone: +49-7141-9730230
Fax: +49-7141-9730240
sales@npielectronic.com
www.npielectronic.com

North America:

ALA Scientific Instruments
Phone: +1-631-393-6401
Fax: +1-631-393-6407
sales@alascience.com
www.alascience.com

Great Britain:

Scientifica Limited
Phone: +44-1825-749933
Fax: +44-1825-749934
info@scientifica.uk.com
www.scientifica.uk.com

Switzerland:

Science Products Trading AG
Phone: +41-43-4880561
Fax: +41-43-4880562
info@science-products.com
www.science-products.ch

