

ALA CPM-2

The CPM-2 is a coating and polishing microforge for micropipettes including an inverted microscope.

Using the CPM-2 kit most inverted microscopes can be converted into a tip processing forge . Thus, the CPM-2 kit can turn any unused inverted microscope into a convenient and useful tool.

Please go to [ALA CPM-2 page](#) for complete information on microforge systems.

npi electronic is distributor for ALA products in Europe except Great Britain.



CPM-2 with Scope

Ordering

Part no.	Description
CPM-2	Coating and polishing microforge kit for inverted microscopes
CPM-2w/scope	Coating and polishing microforge, includes inverted microscope

Download

Brochure:

[CPM-2 brochure](#)